TePe'

Product catalogue





Table of contents

- Sustainability
- Interdental cleaning
- 13 **Toothbrushes**
- 15 **Imprinting**
- Specialty products 20
- Areas of use
- 26 Marketing & display material
- Education 28

Areas of use:



Basic oral care suitable for healthy individuals





Suitable for and/or juniors



Suitable for individuals with braces



Suitable for individuals with implants



individuals
with gingivitis/
periodontitis

Packaging:














From a wooden dental stick to relief for dry mouth

What you are holding in your hand is more than a product catalogue; it's a chronicle of innovative TePe products since 1965, with the latest additions being the toothpaste and products for dry mouth.

Flick through almost 60 years of history of helping people lead healthier lives.

The same dedication, ingenuity and collaborative spirit as Henning Eklund demonstrated

when creating our first product, the wooden dental stick, is present in everything we undertake, initiate and create. All for promoting good oral health.









Plastic coated wire and rounded tip for safety and comfort.



Ergonomic handle for a stable grip.

Plastic sourced from renewable materials.

The plastic cap can be used for hygienic protection or to extend the handle.

TePe Interdental Brushes


According to studies, interdental brushes are the most effective tool for removing dental plague between teeth. We have developed our brushes in collaboration with dental experts to certify the highest functionality, comfort and quality.

The wide TePe Interdental Brush range offers an option for everybody – multiple sizes, two

bristle textures and two handle lengths with a straight or angled brush head.

The brushes comply with the international standard for manual interdental brushes (ISO 16409) which stipulates nine sizes, 0-8. We offer colour-coded brushes in all nine sizes to cover the wide span of interdental space varieties.

Original									
Extra soft									
Angle									
Colour	Pink	Orange	Red	Blue	Yellow	Green	Purple	Grey	Black
ISO size (0-8)	0	1	2	3	4	5	6	7	8
PHD (mm)	0.6	0.7	0.9	1.1	1.3	1.6	2.0	2.6	3.4
Wire size (mm)	0.4	0.45	0.5	0.6	0.7	0.8	1.1	1.3	1.5



The measurement between the dotted lines shows the smallest hole which the brush can pass through, passage hole diameter (PHD), using

TePe® Interdental Brush Original

- The user-friendly handle offers a stable grip which enables cleaning with controlled
- · The bristles cover the entire wire for safe and gentle cleaning.
- The smaller brush sizes (0-3) have a flexible neck for increased durability and access. The larger brush sizes have a thicker, bendable wire.













Smaller sizes: Curve the flexible neck for access between the posterior teeth.

Interdental cleaning.

Larger sizes: Curve the wire for access between the posterior teeth.

Furcations.













TePe® Interdental Brush Original is also available in a mixed pack.

TePe® Interdental Brush Original




^{*}Carton pack is only available in limited markets.

TePe[®] Interdental Brush Extra Soft

- Extra soft bristles for gentle and safe cleaning.
- Six sizes, all with plastic coated wire.
- For patients with sensitive teeth, sore gingiva, dry mouth, or after oral surgery*.



TePe® Interdental Brush Extra Soft



^{*} A dental professional should always be consulted regarding brushing after oral surgery.

TePe Angle™


- · Angled brush head for easier access, especially between the posterior teeth.
- · Plastic coated wire for safe use.
- · Long, flat handle for a stable grip.



TePe Angle™

x25











Colour: Orange Size: Extra Small/Small









Colour: Blue **Size:** Medium/Large









Colour: Purple Size: Extra Large








TePe Mini Flosser™

- Pre-loaded dental floss holder for easy and convenient use in narrow spaces.
- · Strong, thin, teflon free floss.
- Handle made from sugar cane and wood fibres.
- · Bite plane for safe and gentle insertion.











Recommending the correct interdental cleaning device

Individuals with spaces between teeth benefit the most from interdental brushes, which effectively clean the concave and irregular proximal area.



TePe EasyPick™ is a good choice for healthy proximal sites, in the case of gingivitis, and as an on-the-go alternative for everyone.



Floss is recommended for tight, healthy proximal sites, where an interdental brush will not fit.

References: 1. Chapple IL et al. J Clin Periodontol. 2015 Apr;42 Suppl 16:571-6. 2. Ustaoğlu G et al. Clin Oral Investig. 2020 Jun;24(6):2121-2127. 3. Sälzer S et al. Periodontol 2000. 2020 Oct:84(1):35-44.

TePe® Dental Floss

Length: 5m/40m

- The expanding floss adapts itself to the dental space and suits different gaps.
- Slides easily between teeth without shredding due to the vegetable wax/avocado oil coating.
- · Floss made from recycled water bottles.
- Fresh mint flavour.









Gently move the floss along each tooth surface.



TePe® Bridge & Implant Floss

Amount: 30 pcs

- Spongy, flexible floss with stiff ends for easy insertion.
- For daily cleaning of implants, bridges and orthodontic appliances.
- · Intended for one-time use.









Implant-supported bridge.



Practice Box

- Practical storage of TePe's interdental brushes in the dental practice.
- 8 removable trays for IDB Original and Extra Soft, and 4 trays for TePe Angle. Each tray holds 15 brushes.
- A useful information chart with brush sizes comes with the box.

IDB Original /Extra soft

TePe Angle™

IDB Original /Extra soft

TePe Angle™

x1



TePe® Gingival Gel

Size: 20 ml

- Contains chlorhexidine (0.2%) to fight bacteria, and fluoride (0.32% NaF) to strengthen the tooth enamel and protect sensitive root surfaces.
- Non-abrasive, without alcohol and colourants.
- · Fresh mint flavour.















TePe® Travel Case

- Practical travel case for TePe[®] Interdental Brush,
 TePe Mini Flosser™ or TePe EasyPick™.
- · Air holes to allow bristles to dry.
- · For hygienic storage during travel.







TePe Select™ Regular

Size: Regular

Bristles: Extra soft, soft & medium

- · Tapered brush head for easy access.
- · End-rounded bristles for gentle and effective cleaning.
- · User-friendly handle design.
- · Made with renewable materials through the mass balance approach.1







TePe Select™ Compact

Size: Compact

Bristles: Extra soft, soft & medium

- · Small, tapered brush head for easy access.
- · End-rounded bristles for gentle and effective cleaning.
- · Suitable for children and adults who prefer a smaller brush head.
- · Made with renewable materials through the mass balance approach.1







TePe Colour™

Size: Regular Bristles: Soft

- · Tapered brush head for easy access.
- · End-rounded bristles for gentle and effective cleaning.
- · User-friendly handle design.
- · Made with renewable materials through the mass balance approach.1















TePe Colour™ Compact

Size: Compact Bristles: Extra soft

- · Small, tapered brush head for easy access.
- · End-rounded bristles for gentle and effective cleaning.
- · Suitable for children or for adults who prefer a smaller brush head.
- · Made with renewable materials through the mass balance approach.1















TePe Mini™

Size: Mini

Bristles: Extra soft & soft

- · Recommended from the first tooth.
- Extra small, tapered brush head and end-rounded bristles for gentle and effective cleaning.
- Available with or without charming prints created by Swedish Isa Form.
- Made with renewable materials through the mass balance approach.¹













TePe Kids™

Size: Compact
Bristles: Extra soft & soft

- Toothbrush with colourful print, suitable from the age of 3 years.
- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- Made with renewable materials through the mass balance approach.¹













TePe Supreme™

Size: Regular Bristles: Soft

- · Unique two-level bristle design for excellent reach.
- · Easy-grip handle with thumb pad for extra comfort.
- · Made with renewable materials through the mass balance approach.1











TePe Supreme™ Compact

Size: Compact Bristles: Soft

- · Unique two-level bristle design for excellent reach and easy-grip handle with thumb pad.
- · Suitable for children or for adults who prefer a smaller brush head.
- · Made with renewable materials through the mass balance approach.1













Size: Regular

Bristles: Extra soft (pink), soft (yellow) & medium (blue)

- · Long brush tip and end-rounded bristles for effective and gentle cleaning.
- · Easy-grip handle with thumb pad for extra comfort.
- · Made with renewable materials through the mass balance approach.1











Longer tip for improved access to back teeth.







TePe® Extra Grip

- · Provides a comfortable, stable grip which helps patients with reduced manual dexterity.
- · Fits all of our toothbrushes except Supreme, Nova and Choice.
- · Easy to insert and remove by tapping the toothbrush handle on a hard surface.
- · Made with renewable materials through the mass balance approach.1
- · Weighs only 30 grams.









TePe PlagSearch™

Size: 10 tablets/per pack

- · A disclosing product for quick and simple self-check identification of plaque.
- · Two-tone colour indicators.
- · Teaching aid for better oral health.





Differentiates new (red) and old (blue) plaque.



















- · Small, practical travel cap.
- · Air-holes ensure quick drying of the brush.
- · Suits all TePe toothbrushes and special brushes (except the Denture Brush).
- · Made with renewable materials through the mass balance approach.1











TePe Compact Tuft™

- Optimal for precision cleaning.
- Dome-shaped brush head with dense, rounded tuft.
- · Ideal for cleaning along the gum line, on the occlusal surfaces of emerging teeth, or impacted wisdom teeth.
- Made with renewable materials through the mass balance approach.1











TePe Universal Care™

- · A unique angle provides easy access to the lingual and palatal surfaces of teeth and implants.
- · Facilitates cleaning in narrow and hard-toreach areas.
- · End-rounded bristles for gentle cleaning.
- · Made with renewable materials through the mass balance approach.1















TePe® Implant Orthodontic Brush

- · Extra slim brush head for improved access.
- · Optimal for cleaning buccal implant surfaces or orthodontic appliances.
- · Long and thin neck for better access.
- · Made with renewable materials through the mass balance approach.1













TePe Interspace™

Bristles: Soft (yellow) & medium (blue)

- · Ideal for cleaning furcations, implants, fixed braces or behind the last tooth.
- · Brush tip can be attached to either side of the handle.
- · Comes with 12 replaceable brush tips.
- · Made with renewable materials through the mass balance approach.1

















TePe Special Care™

Size: Regular (red) & Compact (blue) Bristles: Ultra soft

- 12,000 (Regular) and 9,000 (Compact) ultra soft bristles for delicate cleaning.
- First brush stage after oral surgery.*
- For patients with extremely sore and delicate tissues, e.g. as a result of trauma, radiotherapy treatment or after oral surgery.
- · Made with renewable materials through the mass balance approach.1





* A dental professional should always be consulted regarding brushing after oral surgery.



TePe Gentle Care™

Size: Regular Bristles: Super soft

- · Super soft filaments for gentle cleaning.
- · Second brush stage after oral surgery.*
- For patients with dry mouth, sensitive root surfaces or sensitive oral tissues.
- · Made with renewable materials through the mass balance approach.1





TePe® Denture Care

Size: Regular

- · Effective and safe cleaning of removable dentures.
- · Brush head can be bent backwards for improved grip and accessibility.
- · Designed to resemble a regular toothbrush.
- · Made with renewable materials through the mass balance approach.1





Long, strong filaments for efficient cleaning of dentures.



TePe® Tongue Cleaner

- · Three blades for a triple cleaning effect.
- Helps prevent bad breath (halitosis).
- Removes bacteria in a few strokes.
- Made with renewable materials through the mass balance approach.1

















TePe® Dental Model

- For patient demonstration and instruction.
- Shows crowding, various interdental spaces, furcation, an erupting wisdom tooth and tooth loss.
- Can be marked with a non-permanent pen for instructional use.



Delivered in a brown hard paper box.



TePe® Implant Care Kit*

A collection of carefully selected maintenance products. Effective plaque removal in and around implants made simple.

Products included:

- TePe Compact Tuft™
- TePe Implant Orthodontic™ Brush
- TePe Universal Care™
- TePe® Bridge & Implant Floss (30 pcs)
- TePe® Interdental Brushes Extra Soft mixed pack
- Patient Information Leaflet





TePe® Orthodontic Kit*

The kit features carefully selected self-care products and readily accessible tips and advice for orthodontic patients.

Products included:

- TePe Supreme™ Compact
- TePe Compact Tuft™
- TePe® Interdental Brush, ISO size 4 and 5
- TePe Angle™ Interdental brush, ISO size 4
- TePe EasyPick™, XS/S and M/L
- TePe® Orthodontic wax
- Patient Information Leaflet







^{*}Currently available in selected markets.

Healthy smiles for all mouths

On this page, you will find recommendations of products for daily oral care for everyone, from the first baby tooth to special needs like oral care for braces and implants.



Paedodontic care

In order to promote a healthy mouth throughout life, we recommend to start brushing when the first tooth emerges. Interdental cleaning is often recommended from the teenage years onwards.





Daily care


The basis for maintaining healthy gums is daily interdental cleaning and toothbrushing. Professionals have a key role in motivating the patient and recommending the proper products.





Implant care


Professional guidance on suitable products and cleaning techniques is important to prevent the development of peri-implant mucositis and peri-implantitis.





Orthodontic care

Fixed orthodontic appliances can lead to a quick build-up of plaque around brackets, bands, and wires. Systematic self-care routines are needed to prevent gingivitis and caries.



Marketing & display material

TePe offers a range of information and demonstration materials to aid oral hygiene advice. Our product racks are perfect for display and storage.



1 Prescription Pad

2 Consumer Leaflets

3 Pharmacy Guide

4 Interdental Brush Guide

TePe[°]

5 Product Catalogue

Item no: 116INT

Caring for your healthy smile: 186INT
Caring for your implants: 172INT
Caring for your braces: 173INT
Caring for your children's teeth: 178INT
Caring for your gums: 184INT
Caring for your mouth throughout life: 1

Item no: 1691N1

Item no: 132IN



Mini Display

Item no: 931 031 Designed for 6 toothbrushes and 9 interdental brushes. Clear acrylic material.



Paper Bag

Item no: 256 INT Height: 32.5 cm. Width: 23 cm. Depth: 10 cm.



Floor Stand, Metal

Item no: 931 054 Floor stand on wheels with 21 hangers, 1 shelf and 2 leaflet holders. Height: 164.5 cm. Width: 50 cm. Depth: 37.5 cm.



Toothbrush Display

Item no: 931 029 Display for 13 toothbrushes (the Select range), Clear acrylic material. Dishwasher proof.



Interdental **Brush Display**

Item no: 931 035 Suitable for 15 interdental brushes. Clear acrylic material.



Bathroom Stand

Item no: 932 048 Suitable for 4 interdental brushes. Clear acrylic material.

EasuPick

Dispenser Box

Item no: 242 291 (XS/S)

242 296 (M/L)

242 271 (XL)

Dispenser box for EasyPick

samples. 100 x 2-pack

included.



Vanity Bag

Item no: 227 Can be filled with products and information for your patient. Partly from renewable material. Height: 14 cm. Width: 25 cm.



Sample Box

Item no: 152 Filled with a variety of product samples and information. Height: 35.5 cm. Width: 23.5 cm. Depth: 5 cm.



EasuPick Displau **Demonstration Case**

Item no: 221INT Cardboard table display suitable for 3 SRPs (shelf ready packaging) with EasyPick. Height: 25 cm. Width: 27 cm. Depth: 16 cm.



Item no: 952 056

Prepacked with the majority of our products to show the variation of the product range. Includes extra space for marketing material and samples. Closes with a zipper. Carry it with ease via the shoulder strap or handle.



Carton Box

Item no: 253 INT Can be filled with products and information for your patient. Height: 23.5 cm. Width: 12 cm. Depth: 2.3 cm.

Medium Table

Item no: 932 033

Table display in cardboard

with 8 hangers. Height: 63.5 cm.

Width: 36.5 cm. Depth: 16 cm.

Display



Mini Flosser Dispenser Box

Item no: 25098 Dispenser box for Mini flosser samples. Height: 16 cm. Width: 16 cm. Depth: 9 cm.



Wall Display, Metal

Item no: 931 053 Wall display with 12 hangers and a shelf. (Plastic containers available separately, item no. 931 038). Height: 80 cm. Width: 50 cm. Depth: 24.5 cm.



Table Display, Metal

Item no: 931 055 Table display with rotating function, 8 hangers. Height: 53 cm. Width: 20 cm. Depth: 34 cm.

Product Presentation Case

Item no: 952 054 Case featuring the major part of the TePe range.





Caring about people's health means that we keep developing solutions for long-term well-being and a better quality of life. The growing scientific support for the link between oral and general health makes creating positive change more relevant than ever. That is why education and awareness is so important.



Articles



Even



ectures.



Sharing materials



Tutorials



Newsletters

We invite you to take part of TePe's knowledge platform. You will find educational events, like webinars, seminars, tutorials and newsletters. Get in-depth knowledge, stay current with the latest research and participate in engaging discussions for your personal development and the benefit of all.

Learn more at **tepe.com**

TePe Munhygienprodukter AB

Bronsåldersgatan 5 | 213 76 Malmö, Sweden Tel +46 (0)40 670 11 00 | Fax +46 (0)40 670 11 11 info@tepe.com | www.tepe.com

All trademarks herein are the property of TePe Munhygienprodukter AB.



TePe is certified according to the international quality standard ISO 9001 and the environmental standard ISO 14000.

